

# User Manual

Instructions for Use

## CRISPincette™ MultiCancer



**Please read all contents in the product manual before use.**

**Research Use Only. Not for use in diagnostic procedures.**

# Table of Contents

1. Product Information.....	4
1.1. Product Description.....	4
1.2. Product Component and Storage Condition .....	4
1.3. Required Equipment and Material (Not Provided) .....	5
2. Testing Protocol.....	6
2.1. Important Notes.....	6
2.2. Time Required for 12 Samples .....	6
2.3. Sample Preparation .....	6
2.4. Library Preparation .....	7
Step 1. 1 <sup>st</sup> PCR.....	7
Step 2. <i>In vitro</i> Cleavage (IVC).....	8
Step 3. Size Selection .....	9
Step 4. 2 <sup>nd</sup> PCR.....	10
Step 5. Size Selection .....	11
Step 6. Library Pooling .....	12
Step 7. Library QC.....	12
2.5. Sequencing (Refer to manufacturer's guidelines) .....	12
Step 8. Denaturation and Dilution Library.....	12
Step 9. Sequencing .....	13
3. Appendix .....	14
A. Index Information .....	14
B. Sample Sheet .....	16
C. Troubleshooting Guide .....	17
D. Precautions for Use .....	18
4. Documentation and Support.....	19

## 1. Product Information



**Important:** Before using this product, read and understand the information in the “D. Precautions for Use” of the **Appendix** in this document.

### 1.1. Product Description

**CRISPincette™ MultiCancer** is a Next-Generation Sequencing (NGS) kit designed to detect mutations in cfDNA from the blood for cancer risk assessment such as lung, breast, colorectal, prostate, gastric, liver, ovarian, and pancreatic cancers. Utilizing GeneCker's ultra-precision CRISPR-Cas9 system, this kit effectively detects even low-frequency mutations by selectively cutting wild-type DNA and amplifying mutant DNA. This product is designed to analyze 23 hotspot mutations across 14 genes, including AKT, APC, BRAF, BRCA1, CTNNB1, EGFR, ERBB2, IDH1, KRAS, MYCN, NRAS, PDGFRA, PIK3CA, and TP53. Results from this kit are intended for clinical reference only. Clinicians should comprehensively interpret the test results in the context of the examinee's condition, drug indications, treatment plan, and other laboratory test indicators.

Major Target List		
AKT1 E17K	APC E1309Dfs*4	BRAF V600E
BRCA1 E23Vfs*17	CTNNB1 T41A	EGFR C797S
EGFR E746_A750del	EGFR L858R	EGFR T790M
ERBB2 V777L	IDH1 R132H/C	KRAS G12D
MYCN P44L	NRAS G12D	PDGFRA D842V
PIK3CA E545K	PIK3CA H1047R	TP53 R175H
TP53 M237I	TP53 R248Q	TP53 R273H
TP53 R337H	TP53 Y220C	

### 1.2. Product Component and Storage Condition



**Important:** Upon arrival, inspect all consumables and contact Warranty Service team if any of the components have been damaged during shipping, store all components under the recommended condition and in an upright position.

Kit	Type	Kit Components	Cap Color	Quantity	Volume (μL)	Storage
<b>Box</b>	tube	● Target primer set1	Green	1	120	-20°C
		● Target primer set2	Yellow	1	120	
		○ 2X PCR mix	White	2	2,400	
		● Remov Enzyme mix set1	Green	1	192	
		● Remov Enzyme mix set2	Yellow	1	192	
		● Remov RXN buffer	Red	1	384	
		● Stabilizer	Red	1	96	
		● 10X STOP buffer	Red	1	96	
<b>Index primer set</b>	Plate	Index primer set (Index 01~96)	-	1	2.5 μL/well	-20°C

## &lt; Index Information &gt;

	1	2	3	4	5	6	7	8	9	10	11	12
A	Index 01	Index 09	Index 17	Index 25	Index 33	Index 41	Index 49	Index 57	Index 65	Index 73	Index 81	Index 89
B	Index 02	Index 10	Index 18	Index 26	Index 34	Index 42	Index 50	Index 58	Index 66	Index 74	Index 82	Index 90
C	Index 03	Index 11	Index 19	Index 27	Index 35	Index 43	Index 51	Index 59	Index 67	Index 75	Index 83	Index 91
D	Index 04	Index 12	Index 20	Index 28	Index 36	Index 44	Index 52	Index 60	Index 68	Index 76	Index 84	Index 92
E	Index 05	Index 13	Index 21	Index 29	Index 37	Index 45	Index 53	Index 61	Index 69	Index 77	Index 85	Index 93
F	Index 06	Index 14	Index 22	Index 30	Index 38	Index 46	Index 54	Index 62	Index 70	Index 78	Index 86	Index 94
G	Index 07	Index 15	Index 23	Index 31	Index 39	Index 47	Index 55	Index 63	Index 71	Index 79	Index 87	Index 95
H	Index 08	Index 16	Index 24	Index 32	Index 40	Index 48	Index 56	Index 64	Index 72	Index 80	Index 88	Index 96

## 1.3. Required Equipment and Material (Not Provided)

Besides the components of the kit, the following equipment and materials are also required.

## 1) List of Equipment and Consumables

	Equipment	Supplier / Catalog No.
1	MiSeq or MiSeq Dx	Illumina
2	96-well Thermal cycler	General laboratory supplier
	Digital Electrophoresis System (Can be selected from the equipment below)	
	Agilent 2100 Electrophoresis Bioanalyzer	Agilent Technologies, Cat # G2939AA
	- Agilent High Sensitivity DNA Kit	Agilent Technologies, Cat # 5067-4626
	Agilent 4150 TapeStation	Agilent Technologies, Cat # G2992AA
3	(or Agilent 4200 TapeStation)	(Cat # G2991BA)
	- Agilent High Sensitivity D1000 ScreenTape & Reagent	Agilent Technologies, Cat # 5067-5584/5067-5585
	- Agilent D1000 ScreenTape & Reagent	Agilent Technologies, Cat # 5067-5582/5067-5583
	- Agilent Genomic DNA ScreenTape & Reagent	Agilent Technologies, Cat # 5067-5365/5067-5366
	- Agilent Cell-free DNA ScreenTape & Reagent	Agilent Technologies, Cat # 5067-5630/5067-5631
	Magnetic Separation Stand (Can be selected from the equipment below)	
	- 16-tube DynaMag™-2 Magnet (for 1.5 mL tube)	Thermo Fisher Scientific, Cat # 12321D
4	- 6-tube Magnetic Separation Rack (for 1.5 mL tube)	New England Biolabs, Cat # S1506S
	- DiaMag02 Magnetic Rack (for 1.5 mL tube)	Diagenode, Cat # B04000001
	- Magnet Stand-96 (DynaMag-96 side) (for 0.2 mL tube)	Thermo Fisher Scientific, Cat # 12331D
	Fluorometer (Can be selected from the equipment below)	
	Qubit® 3.0/4.0 Fluorometer	Thermo Fisher Scientific, Cat # Q33216/Q33226
	- Qubit® Assay tubes	Thermo Fisher Scientific, Cat # Q32856
5	- Qubit® dsDNA Quantitation, High Sensitivity	Thermo Fisher Scientific, Cat # Q32851
	- Qubit® 1x dsDNA HS Assay Kit	Thermo Fisher Scientific, Cat # Q33231
	Quantus Fluorometer	Promega, Cat # E6150
	- QuantiFluor® Single-tube Fluorometer Accessories	Promega, Cat # E4942
	- QuantiFluor® ONE ds DNA System	Promega, Cat # E4871
6	Micro centrifuge	General laboratory supplier
7	Vortex mixer	General laboratory supplier
8	Plate centrifuge	General laboratory supplier
9	Cooling rack	General laboratory supplier

## 2) List of Materials

	Other Materials	Supplier / Catalog No.
1	>99% ethanol	General laboratory supplier
2	microtubes, 1.5 mL	General laboratory supplier

3	8-strip tube	General laboratory supplier
4	DEPC-water	General laboratory supplier
5	SPRIselect, 450 mL	Beckman Coulter, Cat # B23319
6	1.0 N NaOH (molecular biology grade)	General laboratory supplier
7	MiSeq Reagent Kit	Illumina
8	50 mL conical tubes	General laboratory supplier
9	Qubit™ assay tubes	Thermo Fisher, Cat # Q32856
10	Pipette & pipette tip (10, 20, 100, 200 and 1,000 µL)	General laboratory supplier
* Filter tip is recommended		

## 2. Testing Protocol

### 2.1. Important Notes

- 1) All reagents used in this protocol should be performed briefly spin down to collect components clinging under the lid and wall of tubes.
- 2) For frozen reagents, they are fully thawed and mixed using a vortex mixer. Briefly spin down.
- 3) Keep all reagents on ice until ready for use.
- 4) To obtain a high concentration of cfDNA, proceed with the elution volume at the minimum volume guaranteed by the manufacturer.
- 5) If all reagents are repeatedly frozen and thawed, it might cause performance deterioration.
- 6) Ensure magnetic beads are sufficiently vortexed before use.
- 7) 85% ethanol for size selection is prepared and used immediately before every washing process.
- 8) The use of filter tips throughout the entire process is recommended.
- 9) Final library loading concentration of MiSeq is 6 pM, and PhiX control is not required.

### 2.2. Time Required for 48 Samples

Part	Step	Time (hrs)	
Library Preparation	1 <sup>st</sup> PCR	1.5	
	1 <sup>st</sup> PCR & IVC	<i>In vitro</i> Cleavage (IVC)	1
		Size Selection	1
	2 <sup>nd</sup> PCR	2 <sup>nd</sup> PCR	1.5
		Size Selection	0.5
		Library Pooling	0.5
Sequencing	MiSeq or MiSeq DX System	-	

### 2.3. Sample Preparation

\* This product uses DNA that meets the summary of recommended DNA quality and quantity standards below.

#### 1) Recommended tool for sample quantification and qualification

- Accurate quantification of DNA is critical for producing a high-quality NGS library. We recommend using a fluorescence assay such as Qubit® Fluorometer for DNA quantification.
- To accurately determine the size distribution of DNA samples and libraries, we recommend using automated electrophoresis tools such as the TapeStation (Agilent) due to its ability to rapidly evaluate size distribution and quality.

## 2) Summary of recommended DNA quality and quantity

Summary of recommended DNA quality & quantity		
cfDNA	Amount of input DNA (Quantitative values based on Qubit® Fluorometer)	≥ 20 ng
	DNA quality	> 70% cfDNA

## 2.4. Library Preparation

Step 1. 1<sup>st</sup> PCR

<b>Materials (Including Kit)</b>	
-	● Target primer set 1
-	● Target primer set 2
-	○ 2X PCR mix
<b>User preparing materials</b>	
-	1.5 mL microtube
-	8-strip tube
-	DEPC-water
-	DNA sample
-	Thermal cycler

- 1) Perform vortexing and spin down after thawing ● Target primer set 1, ● Target primer set 2, ○ 2X PCR mix in a cooling rack on ice.
- 2) Add the DNA sample and each reagent to the 8-strip tube under the conditions shown in the table below. When handling multiple samples, it is recommended to prepare and use a master mix in a 1.5 mL microtube. Prepare separate master mixes for target primer set 1 and set 2 in that case.

Reagent	Volume per sample (μL)
● Target primer set 1 or ● Target primer set 2	2.5
○ 2X PCR mix	12.5
DNA sample	Variable
DEPC-water	Up to 25
<b>Total Volume</b>	<b>25</b>

\* The recommended amount of cfDNA for each primer set is ≥10 ng.

- 3) Perform vortexing and spin down after mixing each reagent.
- 4) Place the 8-strip tube in the thermal cycler preset under the conditions shown in the table below and run the program.

(Lid temperature: 105°C)

Step Description	Temperature (°C)	Time	Cycles
Pre-denaturation	98	1 min	1
Denaturation	98	10 sec	25
Annealing/Extension	65	40 sec	
Hold	4	∞	1

\* Optional: The recommended ramp rate of thermal cycler is 0.5°C/sec.

## Step 2. *In vitro* Cleavage (IVC)

Kit included materials	
-	● Remov Enzyme mix set 1
-	● Remov Enzyme mix set 2
-	● Remov RXN buffer
-	● Stabilizer
-	● 10X STOP buffer
User preparing materials	
-	Thermal cycler
-	1.5 mL microtube
-	8-strip tube

\* Perform the IVC product from each 1<sup>st</sup> PCR product set 1, set 2 (Step 1.). Store the part of 1<sup>st</sup> PCR product that didn't undergo the IVC procedure on ice or in a refrigerator until use for IVC control in the 2nd PCR (Step 4).

- 1) Perform vortexing and spin down after thawing ● Remov RXN buffer, ● Remov Enzyme mix set 1, ● Remov Enzyme mix set 2, ● Stabilizer in cooling rack on ice.
- 2) Add the 1<sup>st</sup> PCR product (Step 1.) and each reagent to the 8-strip tube under the conditions in the table below. When handling multiple samples, it is recommended to prepare and use a master mix in a 1.5 mL microtube. Prepare separate master mixes for Remov Enzyme mix set 1 and set 2 in that case.

Reagent		Volume per Sample (μL)
IVC	1 <sup>st</sup> PCR Product (set 1 or set 2)	1
	● Remov RXN buffer	4
	● Stabilizer	1
	● Remov Enzyme mix set 1 or	4
	● Remov Enzyme mix set 2	
	<b>Total Volume</b>	<b>10</b>

- 3) Perform vortexing and spin down after mixing each reagent.
- 4) Place the 8-strip tube in the thermal cycler preset under the conditions shown in the table below and run the program.  
(Lid temperature: 60°C)

Step Description	Temperature (°C)	Time
<i>In vitro</i> Cleavage (IVC)	45	30 min

- 5) Collect 10 μL of each set 1 and set 2 products into one well of a new 8-strip Tube. At this time, be careful not to mix the IVC control product and the IVC product after the 30 minutes incubation.
- 6) Add 2 μL of ● 10X STOP buffer to IVC product 20 μL.
- 7) Perform vortexing and spin down.



### Step 3. Size Selection

User preparing materials
<ul style="list-style-type: none"><li>- 8-strip tube</li><li>- DEPC-water</li><li>- &gt;99% ethanol</li><li>- magnetic stand</li><li>- SPRI Beads</li></ul>



\* This step is valid only for IVC products (Step 2.).

- 1) Perform vortexing and spin down sufficiently before using the SPRI beads.
- 2) Mix homogeneously by 10 times of pipetting after adding 40 µL of SPRI beads (Ratio: 2X) to approximately 20 µL to the IVC product (Step 2.).
- 3) Incubate samples for 5 minutes at room temperature.
- 4) Prepare fresh 85% ethanol for washing step.
- 5) Place the sample tube on the magnetic stand until the beads are separated in pellet form (about 2-5 minutes).
- 6) Remove and discard the supernatant using a pipette, without disturbing the beads pellet.
- 7) While keeping the 8-strip tubes on the magnetic stand, add 100 µL of 85% fresh ethanol, incubate for 1 minute at room temperature for beads washing, and then carefully remove and discard the supernatant without disturbing the beads using a pipette.
- 8) Remove all remaining ethanol using a pipette once more and incubate for 1 minute at room temperature for ethanol drying. Be cautious not to over dry the beads pellet.
- 9) Add 10 µL of DEPC-water to each well and mix homogeneously by 10 times of pipetting.
- 10) Place the 8-strip tubes on the magnetic stand until the supernatant is clear and the beads are separated in Pellet form (about 1-2 minutes).
- 11) Transfer 10 µL of the clear supernatant to new 8-strip tubes. At this point, ensure that the elution does not contain beads. If beads are present, place the tube back on the magnetic stand and transfer the elution again.

\* If the used volume of DEPC-water was 10 µL, transfer 8 µL of the supernatant to the 8-strip tubes to prevent the beads from being eluted together.

## Step 4. 2<sup>nd</sup> PCR

<b>Materials (Including Kit)</b>
<ul style="list-style-type: none"> <li>- ○ 2X PCR mix</li> <li>- Index primer set Index 01-96 (plate)</li> </ul>
<b>User preparing materials</b>
<ul style="list-style-type: none"> <li>- 1.5 mL microtube</li> <li>- DEPC-water</li> <li>- Thermal cycler</li> </ul>

\* Perform each 1<sup>st</sup> PCR product set 1, set 2 (Step 1.), IVC Size Selection product (Step 3.).

- 1) Perform vortexing and spin down after thawing ○ 2x PCR in a cooling rack on ice.
- 2) Add the IVC control (Step 1.) and size selection product (Step 3.) and each reagent to the 8-strip tube under the conditions shown in the table below. When handling multiple samples, it is recommended to prepare and use a master mix in a 1.5 mL microtube.

IVC control	
Reagent	Volume per Sample (μL)
Mix of 1 <sup>st</sup> PCR product set 1, set 2 (Step 1.)	<b>0.5 + 0.5</b>
○ 2X PCR mix	12.5
DEPC-water	9
Index primer set Index 01-96 (plate)	2.5
<b>Total Volume</b>	<b>25</b>

IVC	
Reagent	Volume per Sample (μL)
IVC size selection product (Step 3.)	<b>1</b>
○ 2X PCR mix	12.5
DEPC-water	9
Index primer set Index 01-96 (plate)	2.5
<b>Total Volume</b>	<b>25</b>

\* Be cautious of the Index primer plate's direction.

- 3) Perform vortexing and spin down after mixing reagents with each 1<sup>st</sup> PCR product set 1, set 2 (Step 1.), and the IVC size selection product (Step 3.).
- 4) Place the 8-strip tube in the thermal cycler preset under the conditions shown in the table below and run the program.

(Lid temperature: 105°C)

Step Description	Temperature (°C)	Time	Cycles
Pre-denaturation	98	1 min	1
Denaturation	98	10 sec	25
Annealing/Extension	65	40 sec	
Hold	4	∞	1

\* Optional: The recommended ramp rate of thermal cycler is 0.5°C/sec.

## Step 5. Size Selection

### User preparing materials

- 1.5 mL microtube
- DEPC-water
- >99% ethanol
- Magnetic Stand
- SPRI Beads

1) Perform vortexing and spin down sufficiently before using the SPRI beads.

2) Mix 2<sup>nd</sup> PCR product (Step 4.) IVC control and IVC separately.

In this process, the IVC control product and the IVC product should be distinguished. However, size selection is performed by mixing samples of several indexes into one 1.5 mL microtube.

3) After adding SPRI beads (Ratio: 0.8X) to the 2<sup>nd</sup> PCR product (Step 4.) as mentioned in the table below, mix homogeneously by pipetting up and down 10 times.

Number of 2 <sup>nd</sup> PCR Product	Volume of 2 <sup>nd</sup> PCR Product (each)	Size Selection Mixture	
		Volume of 2 <sup>nd</sup> PCR Product	Volume of SPRI Beads
n	10 µL	(n x 10) µL	(n x 8) µL
48 (*Maximum)	10 µL	(48 x 10) µL	(48 x 8) µL

4) Incubate samples for 5 minutes at room temperature.

5) Prepare fresh 85% ethanol for the washing step.

6) Place the sample tube on the magnetic stand until the beads are separated in pellet form (about 2-5 minutes).

7) Remove and discard the supernatant using a Pipette, without disturbing the beads pellet.

8) While keeping the 1.5 mL microtube on the magnetic stand, add a sufficient volume of 85% fresh ethanol to cover the pellets. Incubate for 1 minute at room temperature to wash the beads.

9) Carefully remove and discard the supernatant using a pipette, without disturbing the beads.

10) Remove all remaining ethanol using a pipette once more. Be cautious not to over-dry the beads pellet.

11) After adding DEPC-water to a 1.5 mL microtube as mentioned table below, mix homogeneously by pipetting up and down 10 times.

Number of 2 <sup>nd</sup> PCR Product	* Elution Volume of DEPC-water
8-16	30 µL
16-24	50 µL
25-32	70 µL
33-48	100 µL

12) Incubate samples for 1 minute at room temperature.

13) Place the 1.5 mL microtubes on the magnetic stand until the supernatant is clear and the beads are separated in Pellet form (about 1-2 minutes).

14) Transfer the clear supernatant to new 1.5 mL microtubes. At this point, ensure that the elution does not contain beads. If beads are present, place the tube back on the magnetic stand and transfer the elution again.

\* If the used volume of DEPC-water used was 60 µL, transfer 58 µL of the supernatant to the 1.5mL microtubes to prevent the beads from being eluted together.

## Step 6. Library Pooling

### User preparing materials

- 1.5 mL microtube
- DEPC-water
- Fluorometer

- 1) Qubit® Fluorometer is recommended for the concentration of both the IVC control and IVC product from the size selection product (Step 5.). If the Qubit® Fluorometer is not available, use another type of fluorescent dye-based equipment to measure the dsDNA concentration. When using the Qubit® Fluorometer, utilize Qubit® 1x dsDNA HS Assay Kit.

(Qubit® 1x dsDNA HS working solution : Sample = 198 µL : 2 µL)

- 2) Dilute each sample in a 1.5 mL microtube to a 10 ng/µL concentration using DEPC-water.
- 3) Mix 10 µL of the IVC control and 10 µL of the IVC product in a new 1.5 mL microtube.

## Step 7. Library QC

- 1) Qubit® Fluorometer is recommended for the concentration of the pooled library (Step 6.). If the Qubit® Fluorometer is not available, use another type of fluorescent dye-based equipment to measure the dsDNA concentration. When using the Qubit® Fluorometer, utilize Qubit® 1x dsDNA HS Assay Kit.

(Qubit® 1x dsDNA HS working solution : Sample = 190 µL : 10 µL)

- 2) Convert the measured concentration of the library to molar concentration. (Length of dsDNA: ≈ 220 bp)

\*  $\text{Moles dsDNA (mol)} = \text{mass of dsDNA (ng)} / ((\text{length of dsDNA (bp)} \times 615.96 \text{ g/mol/bp}) + 36.04 \text{ g/mol})$

- 3) Dilute the converted molar concentration to 4 nM using DEPC-Water.

## 2.5. Sequencing (Refer to manufacturer's guidelines)

### Step 8. Denaturation and Dilution Library

\* Thaw the reagent cartridge and pre-chilled HT1 Buffer in a water bath at room temperature. Ensure not to exceed the maximum water line. Store them in the refrigerator until it's time for library denaturation.

\* Dilute 1 N NaOH to 0.2 N NaOH using DEPC-treated Water. It is recommended to use the 0.2 N NaOH within 12 hours after dilution for optimal results.

- 1) Perform vortexing and spin down after adding 0.2 N NaOH 5 µL into 4 nM library 5 µL in a 1.5 mL microtube.
- 2) Incubate samples for 5 minutes at room temperature.
- 3) Add 990 µl prechilled HT1 to the tube containing denatured library. The result is 1 ml of a 20 pM denatured library.
- 4) Perform vortexing and spin down.
- 5) Place the library on ice until the next step.
- 6) Dilute to 6 pM by using the table provided below. (PhiX control is not required.)

6 pM	
20 pM library	180 µL
Prechilled HT1 Buffer	420 µL

- 7) Invert the library diluted to the final concentration about 5 times and spin down.
- 8) Place the final prepared library on Ice until it is loaded onto the MiSeq reagent cartridge.

#### Step 9. Sequencing

- 1) Load 600 µL of the final prepared library onto a MiSeq reagent cartridge.
- 2) Set up a MiSeq Sample Sheet using the Illumina Local Run Manager (LRM) system.  
\* It is recommended to use the sample sheet (.csv) provided by GeneCker
- 3) Start the MiSeq run following the Illumina MiSeq control software instructions.  
\* If you are requesting sequencing from another institution, we recommend that the data throughput per index does not exceed 0.2 Gb.

### Warranty and Liability

GeneCker Co., Ltd. is not responsible for problems caused by using test methods other than the one suggested in this product manual. In the event of a problem with the product, the customer can report the problem within 30 days to GeneCker Warranty Service Team.

Tel: +82-2-3499-4247 / Fax: +82-2-3499-4248 / Web: [www.genecker.com](http://www.genecker.com)

### 3. Appendix

#### A. Index Information

Index Primer	17 Index Sequence	15 Index Sequence	State
Index 01	ATTACTCG	TATAGCCT	A1
Index 02	ATTACTCG	ATAGAGGC	B1
Index 03	ATTACTCG	CCTATCCT	C1
Index 04	ATTACTCG	GGCTCTGA	D1
Index 05	ATTACTCG	AGGCGAAG	E1
Index 06	ATTACTCG	TAATCTTA	F1
Index 07	ATTACTCG	CAGGACGT	G1
Index 08	ATTACTCG	GTA CTGAC	H1
Index 09	TCCGGAGA	TATAGCCT	A2
Index 10	TCCGGAGA	ATAGAGGC	B2
Index 11	TCCGGAGA	CCTATCCT	C2
Index 12	TCCGGAGA	GGCTCTGA	D2
Index 13	TCCGGAGA	AGGCGAAG	E2
Index 14	TCCGGAGA	TAATCTTA	F2
Index 15	TCCGGAGA	CAGGACGT	G2
Index 16	TCCGGAGA	GTA CTGAC	H2
Index 17	CGCTCATT	TATAGCCT	A3
Index 18	CGCTCATT	ATAGAGGC	B3
Index 19	CGCTCATT	CCTATCCT	C3
Index 20	CGCTCATT	GGCTCTGA	D3
Index 21	CGCTCATT	AGGCGAAG	E3
Index 22	CGCTCATT	TAATCTTA	F3
Index 23	CGCTCATT	CAGGACGT	G3
Index 24	CGCTCATT	GTA CTGAC	H3
Index 25	GAGATTCC	TATAGCCT	A4
Index 26	GAGATTCC	ATAGAGGC	B4
Index 27	GAGATTCC	CCTATCCT	C4
Index 28	GAGATTCC	GGCTCTGA	D4
Index 29	GAGATTCC	AGGCGAAG	E4
Index 30	GAGATTCC	TAATCTTA	F4
Index 31	GAGATTCC	CAGGACGT	G4
Index 32	GAGATTCC	GTA CTGAC	H4
Index 33	ATTCAGAA	TATAGCCT	A5
Index 34	ATTCAGAA	ATAGAGGC	B5
Index 35	ATTCAGAA	CCTATCCT	C5
Index 36	ATTCAGAA	GGCTCTGA	D5
Index 37	ATTCAGAA	AGGCGAAG	E5
Index 38	ATTCAGAA	TAATCTTA	F5
Index 39	ATTCAGAA	CAGGACGT	G5
Index 40	ATTCAGAA	GTA CTGAC	H5
Index 41	GAATTCGT	TATAGCCT	A6
Index 42	GAATTCGT	ATAGAGGC	B6
Index 43	GAATTCGT	CCTATCCT	C6
Index 44	GAATTCGT	GGCTCTGA	D6

Index 45	GAATTCGT	AGGCGAAG	E6
Index 46	GAATTCGT	TAATCTTA	F6
Index 47	GAATTCGT	CAGGACGT	G6
Index 48	GAATTCGT	GTA CTGAC	H6
Index 49	CTGAAGCT	TATAGCCT	A7
Index 50	CTGAAGCT	ATAGAGGC	B7
Index 51	CTGAAGCT	CCTATCCT	C7
Index 52	CTGAAGCT	GGCTCTGA	D7
Index 53	CTGAAGCT	AGGCGAAG	E7
Index 54	CTGAAGCT	TAATCTTA	F7
Index 55	CTGAAGCT	CAGGACGT	G7
Index 56	CTGAAGCT	GTA CTGAC	H7
Index 57	TAATGCGC	TATAGCCT	A8
Index 58	TAATGCGC	ATAGAGGC	B8
Index 59	TAATGCGC	CCTATCCT	C8
Index 60	TAATGCGC	GGCTCTGA	D8
Index 61	TAATGCGC	AGGCGAAG	E8
Index 62	TAATGCGC	TAATCTTA	F8
Index 63	TAATGCGC	CAGGACGT	G8
Index 64	TAATGCGC	GTA CTGAC	H8
Index 65	CGGCTATG	TATAGCCT	A9
Index 66	CGGCTATG	ATAGAGGC	B9
Index 67	CGGCTATG	CCTATCCT	C9
Index 68	CGGCTATG	GGCTCTGA	D9
Index 69	CGGCTATG	AGGCGAAG	E9
Index 70	CGGCTATG	TAATCTTA	F9
Index 71	CGGCTATG	CAGGACGT	G9
Index 72	CGGCTATG	GTA CTGAC	H9
Index 73	TCCGCGAA	TATAGCCT	A10
Index 74	TCCGCGAA	ATAGAGGC	B10
Index 75	TCCGCGAA	CCTATCCT	C10
Index 76	TCCGCGAA	GGCTCTGA	D10
Index 77	TCCGCGAA	AGGCGAAG	E10
Index 78	TCCGCGAA	TAATCTTA	F10
Index 79	TCCGCGAA	CAGGACGT	G10
Index 80	TCCGCGAA	GTA CTGAC	H10
Index 81	TCTCGCGC	TATAGCCT	A11
Index 82	TCTCGCGC	ATAGAGGC	B11
Index 83	TCTCGCGC	CCTATCCT	C11
Index 84	TCTCGCGC	GGCTCTGA	D11
Index 85	TCTCGCGC	AGGCGAAG	E11
Index 86	TCTCGCGC	TAATCTTA	F11
Index 87	TCTCGCGC	CAGGACGT	G11
Index 88	TCTCGCGC	GTA CTGAC	H11
Index 89	AGCGATAG	TATAGCCT	A12
Index 90	AGCGATAG	ATAGAGGC	B12
Index 91	AGCGATAG	CCTATCCT	C12
Index 92	AGCGATAG	GGCTCTGA	D12
Index 93	AGCGATAG	AGGCGAAG	E12

Index 94	<b>AGCGATAG</b>	<b>TAATCTTA</b>	F12
Index 95	<b>AGCGATAG</b>	<b>CAGGACGT</b>	G12
Index 96	<b>AGCGATAG</b>	<b>GTACTGAC</b>	H12

## &lt; Index Information &gt;

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	Index 01	Index 09	Index 17	Index 25	Index 33	Index 41	Index 49	Index 57	Index 65	Index 73	Index 81	Index 89
<b>B</b>	Index 02	Index 10	Index 18	Index 26	Index 34	Index 42	Index 50	Index 58	Index 66	Index 74	Index 82	Index 90
<b>C</b>	Index 03	Index 11	Index 19	Index 27	Index 35	Index 43	Index 51	Index 59	Index 67	Index 75	Index 83	Index 91
<b>D</b>	Index 04	Index 12	Index 20	Index 28	Index 36	Index 44	Index 52	Index 60	Index 68	Index 76	Index 84	Index 92
<b>E</b>	Index 05	Index 13	Index 21	Index 29	Index 37	Index 45	Index 53	Index 61	Index 69	Index 77	Index 85	Index 93
<b>F</b>	Index 06	Index 14	Index 22	Index 30	Index 38	Index 46	Index 54	Index 62	Index 70	Index 78	Index 86	Index 94
<b>G</b>	Index 07	Index 15	Index 23	Index 31	Index 39	Index 47	Index 55	Index 63	Index 71	Index 79	Index 87	Index 95
<b>H</b>	Index 08	Index 16	Index 24	Index 32	Index 40	Index 48	Index 56	Index 64	Index 72	Index 80	Index 88	Index 96

**B. Sample Sheet**

Use the sample sheet (.csv) **provided by GeneCker** or prepare a proper sample sheet for Illumina Local Run Manager (LRM) by following the Illumina Experiment Manager software guide.

\* Rules for writing sample sheet: Control - Sample Name\_Visit number\_Control / IVC - Sample Name\_Visit number\_IVC

**Example.**

1	[Header]						
2	Experiment Name	Test1	①				
3	Date	0000-00-00	②				
4	Module	GenerateFASTQ - 2.0.0					
5	Workflow	GenerateFASTQ					
6	Library Prep Kit	Custom	③				
7	Chemistry	Amplicon					
8	[Reads]						
9		151					
10		151					
11	[Data]						
12	Sample_ID ④	Description	I7_Index_ID	index	I5_Index_ID	index2	Sample_Project
13		1	A701	ATTACTCG	A501	tatagcct	
14		2	A701	ATTACTCG	A502	atagaggc	
15		3	A701	ATTACTCG	A503	cctatcct	
16		4	A701	ATTACTCG	A504	ggctctga	
17		5	A701	ATTACTCG	A505	aggcgaag	
18		6	A701	ATTACTCG	A506	taatctta	
19		7	A701	ATTACTCG	A507	caggacgt	
20		8	A701	ATTACTCG	A508	qtactgac	

① Insert experiment title.

② Insert NGS run date.

③ Insert 'Custom'. (\* The sample sheet provided by GeneCker has 'Custom' entered.)

④ Insert the Sample\_ID as the example below.

\* Patients' ID\_Tissue or Plasma\_Visit Number\_Control or IVC

\*\* ' ' or ' - ' cannot be enter in sample sheet.

Ex) Patient1234\_Plasma\_V1\_Control (If used index 4, enter 4 instead)

Patient1234\_Plasma\_V1\_IVC (If used index 8, enter 8 instead)

☛ If necessary, the "Sample Sheet" form can be sent by e-mail ([info@genecker.com](mailto:info@genecker.com)) to request.



## C. Troubleshooting Guide

Troubles	Probable Cause	Resolution
Insufficient library volume	Not following the manufacturer's recommended PCR conditions may result in improper amplification of amplicons, leading to insufficient amounts of the library being produced in subsequent steps.	Adhere to the manufacturer's recommended PCR conditions.
	Improper temperature of the PCR lid means DNA evaporation may occur, potentially resulting in insufficient production of the library.	
	If there is an issue with the temperature sensor of the PCR equipment, the proper maintenance of temperature within the PCR device may be compromised, potentially leading to a failure in the correct amplification process.	Regularly inspect the equipment to ensure its proper functioning.
In case the VAF (Variant Allele Frequency) doesn't change value is low or not after the <i>In vitro</i> Cleavage reaction	<b>Short Reaction Time:</b> If the reaction time for <i>in vitro</i> Cleavage is too short, there may be limited cleavage of wild-type DNA, resulting in difficulty obtaining enough of the desired target genes	Follow the manufacturer's recommended <i>in vitro</i> Cleavage conditions and consider re-reaction if necessary.
	<b>Excessive Reaction Time:</b> If the reaction time for <i>in vitro</i> Cleavage exceeds the manufacturer's recommended duration, excessive cleavage of the amplicon may occur.	
	<b>Delayed Stop reaction:</b> If the addition of ●10X STOP buffer is delayed after the <i>in vitro</i> Cleavage reaction, excessive cleavage of the amplicon may occur.	It is recommended to add the ●10X STOP buffer within 5 minutes after the <i>in vitro</i> Cleavage reaction.
	Using the Remov Enzyme mix with an inappropriate set configuration may lead to excessive cleavage of the amplicon.	Perform <i>in vitro</i> cleavage using a Remov Enzyme mix set with the same as the primer set.
If there is a loss of beads during the Size Selection process	Beads may be sucked into the pipette tip during the removal of supernatant.	Remove the supernatant after the magnetic beads are completely bound to a magnetic stand.
In case of inaccurate volume measurement of the beads	Inadequate vortexing of beads may lead to the capture of a non-specific size library due to inaccurate ratios during the process.	Ensure thorough vortexing of the beads and prevent any beads from remaining on the inner and outer walls of the pipette tip.
If beads do not fully dissolve in DEPC-Water after an ethanol wash	If the beads were over-dried.	Pipette until the beads are completely dissolved.
		Immediately dissolve the beads after an ethanol wash to prevent over drying.

#### D. Precautions for Use

1. This product may not be used for purposes other than research.
2. This product is disposable and cannot be reused.
3. It is recommended to read the following precautions before using this product.
4. This product is designed to use Next-Generation Sequencing (NGS) and requires the use of MiSeq or MiSeq Dx equipment, which are base sequencing analyzers.
5. Prepare a master mix of all reagents on ice.
6. Since the NGS library preparation step is sensitive to pipetting errors, it is recommended to use a calibrated pipette within 1 year.
7. It is recommended to use filter tips to prevent cross-contamination.
8. Since the quality of the DNA sample affects the results, it is recommended to use high-quality DNA. Severely damaged DNA can cause the failure of library construction.
9. The concentration of DNA should be measured using the fluorescent dye method.
10. Contamination by DNA and RNA other than the sample may affect the quantification of DNA extracted from the sample, so keep the surrounding environment clean to avoid cross-contamination.
11. Contamination of test reagents and equipment, reaction temperatures, or storage conditions different from those recommended may affect test results.
12. If you need help with GeneCker's product, please contact the GeneCker Co., Ltd. Support at **info@genecker.com**.












## 4. Documentation and Support

### DISCLAIMER

TO THE EXTENT ALLOWED BY LAW, **GeneCker Co., Ltd.** WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

### Symbol Descriptions

The following table describes the symbols on the shipment packaging, consumable, or consumable packaging.

Symbol	Description	Symbol	Description
	Indicates the manufacturer.		Storage temperature ranges in degrees Celsius. Store the consumables within the indicated range.
	Date of manufacture.		Caution.
	Indicates the part number so that the consumable can be identified.		The date the consumable expires. For best results, use the consumable before this date.
	Indicates the batch code to identify the manufacturing batch or a lot of the consumables.		Contains sufficient for <n> tests.
	Consult instructions for use.		Do not reuse.
	The intended use is Research Use Only (RUO).		

### Revision History

Document Number	Revision	The Date of Revision	Contents of Revision
GC-IFU-N102R-EN	0	01 July 2024	Initial release

## Customer and Technical support :)

For Customer and technical assistance, contact GeneCker Support.

[info@genecker.com](mailto:info@genecker.com)

[www.genecker.com](http://www.genecker.com)



**GeneCker Co., Ltd.**

**#1202~1204, HausD SeJong Tower, 26, Seongsuil-ro 10-gil, Seongdong-gu, Seoul, 04793 REPUBLIC OF KOREA**

**Tel: 02-3499-4247 / Fax: 02-3499-4248 / Web: [www.genecker.com](http://www.genecker.com)**

GeneCker Co., Ltd. is responsible for the production of products related to the operations mentioned in this user manual.  
The information in this user manual is subject to change without notice.