

CRISPincette™ MultiCancer Quick Guide

Preparation for Assay

- Before starting, spin down all reagents briefly to collect components clinging under the lid and keep them on ice until ready for use.
- In the case of frozen reagents, they are fully thawed and mixed using a vortex mixer. Briefly spin down and keep on ice until ready for use.
- SPRI beads must be sufficiently vortexed before use. It should appear homogenous and consistent in color.
- The 85% ethanol used for Size Selection is prepared and used immediately before the washing process.
- Using a filter tip during all experimental processes is recommended.

Step 1. 1st PCR (Total volume : 25 μ L, each)

1) Add DNA sample and each reagent in 8-strip tube under the conditions shown in the table below.

*** Perform PCR by adding the Target primer sets to each separated tube.**

Reagent	Volume per Sample
● Target primer set 1	2.5 μ L to each tube
● Target primer set 2	
○ 2X PCR mix	12.5 μ L
DNA sample	Variable (Maximum 10 μ L)
DEPC-Water	Up to 25 μ L
Total Volume 25 μ L	

2) Vortex and spin down.

3) Place the 8-strip tube in the pre-set thermal cycler under the conditions shown in the table below and run the program.
(Lid temperature: 105 °C and ramp rate 0.5 °C/sec)

Step description	Temperature	Time	Cycles
Pre-denaturation	98 °C	1 min	1
Denaturation	98 °C	10 sec	25
Annealing/Extension	65 °C	40 sec	
Hold	4 °C	∞	1

Step 2. IVC (Total volume : 10 μ L, each)

1) Add 1st PCR product (Step 1.) and each reagent in 8-strip tube under the conditions shown in the table below.

*** Separately perform IVC on set 1 and set 2 of 1st PCR products (Step 1.) by adding the Remov Enzyme mix sets to the corresponding 1st PCR product sets. The remaining 1st PCR products are used as 'Control (Non-IVC)' in Step 4.**


Reagent	Volume per Sample
1 st PCR product (set 1 or set 2)	1 μ L to each tube
● Remov RXN buffer	4 μ L
● Stabilizer	1 μ L
● Remov Enzyme mix set 1	4 μ L to each tube
● Remov Enzyme mix set 2	
Total Volume 10 μ L	

2) Vortex and spin down.

3) Place the 8-strip tube in the pre-set thermal cycler under the conditions shown in the table below and run the program.
(Lid temperature: 60 °C)

Step description	Temperature	Time
In vitro Cleavage (IVC)	45 °C	30 min


4) After the 30 minutes reaction, collect 10 μ L of each IVC product set 1 and set 2 into one well of a new 8-strip tube. And add 2 μ L ● 10X STOP buffer to each 8-strip tube.



Step 3. Size Selection (with a ratio of 2X)

*** Only perform 'Size Selection' on the IVC products (Step 2.).**

1) Add 40 μ L of SPRI beads (2X ratio) to approximately 20 μ L of IVC product (Step 2.). Then, perform 'Size Selection' according to the following steps.



CRISPincette™ MultiCancer Quick Guide

Step 4. 2nd PCR (Total volume : 25 μ L, each)

* Perform PCR on the Control (1st PCR product 1 and 2, Step 1.) and on the IVC Size Selection product (Step 3.), respectively.

1) Add Control and IVC Size Selection product and each reagent in 8-strip tube under the conditions shown in the table below.

Reagent	Volume per Sample
Control (Mix 1 st PCR product 1 and 2)	0.5 + 0.5 μ L
○ 2X PCR mix	12.5 μ L
DEPC water	9 μ L
Index primer set index 01-24 (plate)	2.5 μ L
Total Volume 25 μ L	

Reagent	Volume per Sample
IVC Size Selection product	1 μ L
○ 2X PCR mix	12.5 μ L
DEPC water	9 μ L
Index primer set index 01-24 (plate)	2.5 μ L
Total Volume 25 μ L	

2) Vortex and spin down.


3) Place the 8-strip tube in the pre-set thermal cycler under the conditions shown in the table below and run the program.
(Lid temperature: 105 °C and ramp rate 0.5 °C/sec)

Step description	Temperature	Time	Cycles
Pre-denaturation	98 °C	1 min	1
Denaturation	98 °C	10 sec	25
Annealing/Extension	65 °C	40 sec	
Hold	4 °C	∞	1

Step 5. Size Selection (with a ratio of 0.8X)

1) Separately mix the 2nd PCR products (Step 4.) of the Controls and IVC products into each 1.5 mL microtube. Add SPRI beads with a ratio of 0.8X to the 1.5 mL microtube under the conditions shown in the table below.

Size Selection mixture			
Number of 2 nd PCR product	Volume of 2 nd PCR product (each)	Volume of 2 nd PCR product	Volume of SPRI beads
n	10 μ L	(n \times 10) μ L	(n \times 8) μ L
12 (Maximum)	10 μ L	(12 \times 10) μ L	(12 \times 8) μ L



* Final elution volume of DEPC-water

Number of 2 nd PCR product	Elution Volume of DEPC-Water
<10	30 μ L
10-12	50 μ L

Step 6. Library Pooling

1) Qubit® Fluorometer equipment is recommended to measure the concentration (ng/ μ L) of Size Selection products (Step 5.).

2) Dilute Size Selection products in new 1.5 mL microtubes to 10 ng/ μ L using DEPC-water.

3) Mix each 10 μ L of the diluted Size Selection products in a new tube.

Step 7. Library QC

1) Qubit® Fluorometer equipment is recommended to measure the concentration (ng/ μ L) of the pooled library.

2) Convert the measured concentration of the pooled library to molar concentration.
(Length of dsDNA : \approx 220 bp)

3) Dilute the pooled library to 4 nM with DEPC-water.

Step 8. Sequencing

1) Please refer to the manufacturer's instructions and recommendations.