

## CRISPincette™ LUNG Quick Guide

### Preparation for Assay

- Before starting, spin down all reagents briefly to collect components clinging under the lid and keep them on ice until ready for use.
- In the case of frozen reagents, they are fully thawed and mixed using a vortex mixer. Briefly spin down and keep on ice until ready for use.
- SPRI beads must be sufficiently vortexed before use. It should appear homogenous and consistent in color.
- The 85% ethanol used for Size Selection is prepared and used immediately before the washing process.
- Using a filter tip during all experimental processes is recommended.

### Step 1. 1<sup>st</sup> PCR (Total volume : 25 $\mu$ L, each)

1) Add DNA sample and each reagent in 8-strip tube under the conditions shown in the table below.

**\* Perform PCR by adding the Target primer sets to each separated tube.**

Reagent	Volume per Sample
● Target primer set 1	2.5 $\mu$ L to each tube
● Target primer set 2	
○ 2X PCR mix	12.5 $\mu$ L
DNA sample	Variable (Maximum 10 $\mu$ L)
DEPC-Water	Up to 25 $\mu$ L
<b>Total Volume 25 <math>\mu</math>L</b>	

2) Vortex and spin down.

3) Place the 8-strip tube in the pre-set thermal cycler under the conditions shown in the table below and run the program.  
(Lid temperature: 105 °C and ramp rate 0.5 °C/sec)

Step description	Temperature	Time	Cycles
Pre-denaturation	98 °C	1 min	1
Denaturation	98 °C	10 sec	30
Annealing/Extension	65 °C	40 sec	
Hold	4 °C	$\infty$	1

### Step 2. IVC (Total volume : 10 $\mu$ L, each)

1) Add 1<sup>st</sup> PCR product (Step 1.) and each reagent in 8-strip tube under the conditions shown in the table below.

**\* Separately perform IVC on set 1 and set 2 of 1<sup>st</sup> PCR products (Step 1.) by adding the Remov Enzyme mix sets to the corresponding 1<sup>st</sup> PCR product sets. The remaining 1<sup>st</sup> PCR products are used as 'Control (Non-IVC)' in Step 4.**


Reagent	Volume per Sample
1 <sup>st</sup> PCR product (set 1 or set 2)	1 µL to each tube
● Remov RXN buffer	4 µL
● Stabilizer	1 µL
● Remov Enzyme mix set 1	4 µL to each tube
● Remov Enzyme mix set 2	
Total Volume 10 µL	

2) Vortex and spin down.

3) Place the 8-strip tube in the pre-set thermal cycler under the conditions shown in the table below and run the program.  
(Lid temperature: 60 °C)

Step description	Temperature	Time
In vitro Cleavage (IVC)	45 °C	30 min


4) After the 30 minutes reaction, collect 10  $\mu$ L of each IVC product set 1 and set 2 into one well of a new 8-strip tube. And add 2  $\mu$ L ● 10X STOP buffer to each 8-strip tube.



### Step 3. Size Selection (with a ratio of 2X)

**\* Only perform 'Size Selection' on the IVC products (Step 2.).**

1) Add 40  $\mu$ L of SPRI beads (2X ratio) to approximately 20  $\mu$ L of IVC product (Step 2.). Then, perform 'Size Selection' according to the following steps.



## CRISPincette™ LUNG Quick Guide

### Step 4. 2<sup>nd</sup> PCR (Total volume : 25 µL, each)

\* Perform PCR on the Control (1<sup>st</sup> PCR product 1 and 2, Step 1.) and on the IVC Size Selection product (Step 3.), respectively.

1) Add Control and IVC Size Selection product and each reagent in 8-strip tube under the conditions shown in the table below.

Reagent	Volume per Sample
Control (Mix 1 <sup>st</sup> PCR product 1 and 2)	0.5 + 0.5 µL
○ 2X PCR mix	12.5 µL
DEPC water	9 µL
Index primer set index 01-48 (plate)	2.5 µL
Total Volume 25 µL	

Reagent	Volume per Sample
IVC Size Selection product	1 µL
○ 2X PCR mix	12.5 µL
DEPC water	9 µL
Index primer set index 01-48 (plate)	2.5 µL
Total Volume 25 µL	

2) Vortex and spin down.


3) Place the 8-strip tube in the pre-set thermal cycler under the conditions shown in the table below and run the program.  
(Lid temperature: 105 °C and ramp rate 0.5 °C/sec)

Step description	Temperature	Time	Cycles
Pre-denaturation	98 °C	1 min	1
Denaturation	98 °C	10 sec	25
Annealing/Extension	65 °C	40 sec	
Hold	4 °C	∞	1

### Step 5. Size Selection (with a ratio of 0.8X)

1) Separately mix the 2<sup>nd</sup> PCR products (Step 4.) of the Controls and IVC products into each 1.5 mL microtube. Add SPRI beads with a ratio of 0.8X to the 1.5 mL microtube under the conditions shown in the table below.

Size Selection mixture			
Number of 2 <sup>nd</sup> PCR product	Volume of 2 <sup>nd</sup> PCR product (each)	Volume of 2 <sup>nd</sup> PCR product	Volume of SPRI beads
n	10 µL	(n × 10) µL	(n × 8) µL
12	10 µL	(12 × 10) µL	(12 × 8) µL
24 (Maximum)	10 µL	(24 × 10) µL	(24 × 8) µL



\* Final elution volume of DEPC-water

Number of 2 <sup>nd</sup> PCR product	Elution Volume of DEPC-Water
<10	30 µL
10-15	50 µL
15-20	70 µL
20-24	100 µL

### Step 6. Library Pooling

1) Qubit® Fluorometer equipment is recommended to measure the concentration (ng/µL) of Size Selection products (Step 5.).

2) Dilute Size Selection products in new 1.5 mL microtubes to 10 ng/µL using DEPC-water.

3) Mix each 10 µL of the diluted Size Selection products in a new tube.

### Step 7. Library QC

1) Qubit® Fluorometer equipment is recommended to measure the concentration (ng/µL) of the pooled library.

2) Convert the measured concentration of the pooled library to molar concentration.  
(Length of dsDNA : ≈ 220 bp)

3) Dilute the pooled library to 4 nM with DEPC-water.

### Step 8. Sequencing

1) Please refer to the manufacturer's instructions and recommendations.